



# Introducing: The Squamous Spectrum

Unsatisfactory and benign/reactive  
samples

Margaret Sage  
March 2021

# Bethesda 2001

Atypical Squamous Cells (ASC)

- of undetermined significance (ASC-US)

**LSIL:** Low-grade Squamous Intraepithelial Lesion

Atypical Squamous Cells (ASC)


- cannot exclude HSIL (ASC-H)

**HSIL:** High-grade Squamous Intraepithelial Lesion

(+/- with features suspicious for invasion)

**Squamous Cell Carcinoma**

## CYTOLOGY



# Squamous abnormalities

Cytology

NZ samples \*

ASC-US

1.8%

LSIL

3.5%

*Low-grade*

ASC-H

0.5%

*High-grade*

HSIL



0.8%

- clinical outcome justifies cytology result categories
- ASC-US + LSIL cytology identifies almost as many HSIL histology cases as ASC-H + HSIL cytology

\*% of satisfactory samples Jan-June 2017 NCSP Monitoring Report 47



Condensation of peripheral cytoplasm

# LSIL: classic koilocytes




Multinucleation

# LSIL: Koilocytes




LSIL: not koilocytes



*ThinPrep*

# ASC-US

## Atypical Squamous Cells of Undetermined Significance




*SurePath*

# HSIL


## High N:C ratio

Hyperchromatic granular chromatin  
Irregular nuclear outline



## Nuclear Variation


Size, shape, chromatin, outline



*ThinPrep*

## Crowded sheets


Thick, disorganised, high N:C  
cells, nuclear variability



*SurePath*




ASC-H  
Atypical Squamous Cells, possible High-grade



*SurePath*



Two suspicious cells... Follow-up = CIN3








*ThinPrep*

Suspicious, some reactive features  
Follow-up = inflammation, cervicitis

# Squamous cell carcinoma




# Squamous cell carcinoma



*ThinPrep*

# REACTIVE Sheet



*ThinPrep*

# UNSATISFACTORY SAMPLES

- A. Rejected specimens – LBC vial leaking, unlabeled.
- B. Specimen examined but unsatisfactory for evaluation

The specimen is **unsatisfactory** for evaluation because....

of **insufficient squamous cells**.

of **poor fixation/preservation**.

**foreign material** obscures the cells.

**inflammation** obscures the cells.

**blood** obscures the cells.

of **cytolysis/autolysis**.

**Bethesda Criteria for adequate cellularity:**  
**Squamous cells well-visualised and well-preserved**

Liquid-based samples: **at least 5,000 cells**

- a minimum of 10 fields counted randomly along a diameter that includes the centre of the preparation

- minimum numbers of cells needed:

**SurePath: 9 cells** per 40X in each of 10 fields

**ThinPrep: 4 cells** per 40X in each of 10 fields

Ref: The Bethesda System for Reporting Cervical Cytology  
3rd Edition. Nayar and Wilbur (eds). Springer.

# Inadequate rates UK

Health Technology Assessment 2015 (Kitchener et al)

<http://www.journalslibrary.nihr.ac.uk/hta/volume-19/issue22#abstract>.

Conclusion: SurePath slides: 15,000 minimum

ThinPrep slides: 5,000 minimum

2013-14 Inadequate rate across UK = 2.4%

Reference:


ABC3 and LBC – Adequate or not?

Duval E. (editorial) *Cytopathology* 2013, **24**, 211-5

# Unsatisfactory rates in New Zealand


Total for NZ samples **July - Dec 2015: 1.3%**

Target: **1-5%** of all LBC samples reported as unsatisfactory







The specimen is **unsatisfactory** for evaluation because of **insufficient squamous cells**.





*Few or no cells*



*Mainly endocervical cells*

The specimen is **unsatisfactory** for evaluation because **blood obscures the cells.**



# Organisms

There are organisms consistent with *Trichomonas vaginalis*

There are fungal organisms morphologically consistent with  
*Candida* species

There is a shift in microbiological flora suggestive of  
bacterial vaginosis


There are bacteria morphologically consistent with  
*Actinomyces* species

There are cellular changes consistent with *Herpes simplex*  
*virus*


# *Trichomonas vaginalis*



Single trichomonads  
insert: flagellum




Large cluster of  
trichomonads




*ThinPrep*

Fungal organisms consistent with *Candida* species






*ThinPrep*




*SurePath*

Shift in bacterial flora suggesting bacterial vaginosis

## Bacteria consistent with *Actinomyces*





*ThinPrep*



*SurePath*

Molding, multinucleation and  
margination of chromatin

Intranuclear viral inclusions

Cell changes consistent with *Herpes simplex* virus



## Reactive/non-neoplastic changes


There are **reactive** cellular changes present.

There are **endometrial cells present in a woman over the age of 40 years.\***

There are **atrophic** cellular changes present.

*\*The presence of endometrial cells in a woman over the age of 40 years can be a normal finding, or seen in association with hormone replacement therapy , or rarely, associated with endometrial pathology including hyperplasia or neoplasia. Please correlate this finding with any symptomatology of uterine pathology, for example abnormal uterine bleeding and refer/investigate appropriately.*





# Reactive squamous cells





# Reactive endocervical cells




# Radiation change




# IUCD cells



# Lymphocytic cervicitis




Loose cluster of lymphocytes




Tingible body macrophages

# Parakeratosis




Squamous pearl



Parakeratotic group showing  
smaller darker nuclei



# Cervical endometriosis



Endometrial cells, histiocytes and blood




# Tubal metaplasia




Columnar endocervical cells with cilia and terminal bars

## Transitional metaplasia



Characteristic longitudinal nuclear grooves

# ATROPHY



See as many cases and images as you can